
MooDLIGHT

TENS

Model:SM9126

ACTIVATE YOUR 24 MONTH WARRANTY

PLEASE VISIT NOW

→<https://moodlightlife.com/warranty>

CUSTOMER SERVICE

support@moodlightlife.com

Mood Light

Please Download User Manual on Below Link:

<https://moodlightlife.com/downloads>

Included in this standard package

- Device x 1pc
- Large pad x 2pcs
- Medium pad x 6pcs
- Small pad x 2pcs
- Electric wire x 3pcs


When use 4 line cable, make sure to use at least 1pc white(positive)+1pc blue line (negative)at the same time.

- USB cable x1pc
- Instruction manual x 1pc
- Bagx1pc

Features

- 2 output with 24 modes
- 20 level intensity
- 20 minutes default timer, can be adjustable from 10-60 minutes
- Rechargeable lithium battery
- Battery indicator to know when need to recharge battery
- AC adapter 110-240v

Illustrations of the device and accessories


1. On/Off switch
2. USB port
3. Battery indicator
4. Time
5. Mode selector
6. A Intensity adjustment (+)
7. Mode adjustment
8. A Intensity adjustment (-)
9. Time adjustment
10. B Intensity adjustment (+)
11. Mode adjustment
12. B Intensity adjustment (-)
13. Pause
14. A output connector
15. B output connector
16. Pad wires
17. Massage pad
18. USB cable

Modes and Function

There are 24 modes to choose from. The “M” button cycles through the 6 main modes, all of these 6 main modes contain different built-in 4 modes.

Main Mode 1:

1) Lateral palm kneading(TENS): It's specially effective for joint pain, sciatica, muscle spasms, and pain in the back, neck, and shoulders.

2) Flat palm kneading(TENS): Enhance function, aid in the healing process, decrease muscle reflex activity, promote relaxation, recreational and well-being.

3) Crosscut fist kneading(EMS): It is a type of alternative medicine consisting of finger and palm pressure, stretches, and other massage techniques.

4) Vertical fist kneading(TENS): It is the manipulation of superficial and deeper layers of muscle and connective tissue using various techniques.

Main Mode 2:

1) Acupuncture(EMS): This mode is used for the points on the Chinese meridian pressure point chart to assist headaches, a variety of ailments as well as muscle and joint pain. Particularly good for fast relief of acute pain.

2) Moxibustion(EMS): Helps relieve pain, relax the mind and muscle, promote sleep.

3) Acupuncture kneading(EMS): Helps muscle reduction, relaxation of muscle spasms, increase range of motion, injury recovery and post surgery recovery.

4) Moxibustion kneading(EMS): It helps chronic conditions such as insomnia Constipation, headaches, and the tension associated with stress.

Main Mode 3:

1) Rubbing(TENS): Provides the user with a scratching and kneading sensation. Releases

unhealthy elements from injured areas and stimulates blood flow and healing.

2) Rubbing reduction(TENS): Is effective in acute and chronic internal organ disorders. Help relieve pain, stiffness, fever, chill, cough, wheeze, nausea and vomiting etc.

3) Tapping(TENS): produce a hard thumping sensation that massage your muscles, which can give you an invigorated and revitalized feeling.

4) Fibrillation massage(TENS): It helps muscle reduction, relaxation of muscle spasms, increase range of motion, injury recovery and post surgery recovery.

Main Mode 4:

1) Lifting-pinching manipulation

(TENS+EMS): It strengthen and firm muscles at abs, back, thighs and arms etc.

2) Vibrating manipulation(TENS+EMS):

For muscle strains, herniated discs, stenosis, spondylosis, spondylolisthesis, and general arthritis pain.

3) High frequency manipulation(TENS):

Auto combination of several modes for relaxation. Full body relaxation after workout: exercises & fitness or after a day of hard work.

4) Fibrillation manipulation(TENS): Help for muscle knots and muscle tightness caused by strain, overworking of your muscles.

Main Mode 5:

1) Shiatsu-Deep(TENS): Help for foot pain caused by plantar fasciitis, scarring, collagen breakdown, and micro-tears, Heel spurs, foot muscles and ligament strains.

2) Thai Massage(TENS): Provide relief from pain and stiffness associated with sports injuries, whiplash, osteoarthritis, muscle spasms, fibromyalgia and strain injuries such as carpal tunnel syndrome.

3) Shiatus(TENS+EMS): Shiatsu can help facilitate healing from many minor ailments, such as headaches and back pain to digestive, sinus and skin problems.

4) Strength Alternate Massage

(TENS+EMS): It can be used to relieve neck pain, shoulder pain.

Main Mode 6:

1) Twinkling Cupping(TENS): Use low frequency vibrations to promote circulation of blood and relieve pain. Specifically use for lower back.

2) Cupping(TENS): For relieve back and neck pains, stiff muscles, anxiety, fatigue, migraines, rheumatism, and even cellulite etc.

3) Fibrillation Cupping(TENS): For muscle strains, herniated discs, stenosis, spondylosis, spondylolisthesis, and general arthritis pain.

4) Cup retaining(TENS): It can be Used for achilles bursitis, sprains, strains, heel pain, bone spurs, tendonitis, ankle impingement, and general arthritis.

Operating instructions

1. Connect the pads to the wire by snapping them on, then remove the protective film.

Note: adding a few drops of water on the pads will increase their longevity.

2. Place pads where the massage sensations

are desired. And plug into one of the ports on the device. For additional relief, use the second port and repeat steps 1&2 for two additional pads.

Make sure both pads on the skin and not overlapping. Make sure your skin is free from any dirt, oil or lotion.

TIPS: When you use two pads (one channel) or four pads(two channels) at the same time, make sure you have placed all the pads on your body, otherwise it won't work.

3. Turn the device on by pressing the on/ off switch on the top of the unit. The unit will be now activated and start in Mode 1.

4. There are six icons on the screen, each icon contains different built-in 4 modes. Press "<" and ">" button to change the six icons (6 main modes). Press the "M" to try 4 different modes within that selection.

Please increase the intensity because it will return to lowest automatically when you change to a new mode.

5. Once you select a mode, gradually increase the A/B intensity by pressing the "+" button on the right/ left side; decrease it by pressing "-" button on the right/ left side of the tens.

Tips: When you plug the wire to A port, you need to press the +/- button on the left side to adjust the intensity.

6. Adjust the tens using time by pressing the T button on the side of the tens. The timer can be increased by 10 minute intervals and the tens will automatically turn off after the set time has finished.

7. Press the round button to pause when change the massage point.

8. Turn the device off by pressing the on/ off switch on the top of the unit If you need to turn off the device during massaging.

9. When battery is low, recharge the battery with 5V-0.5A USB charger 40 minutes.

Precautions

1. Do not apply the massage pads near the heart, on the head,in the pubic region, over Scarred areas, on the throat or over the mouth.

2. Avoid touching the pads when the unit is on.

3. Pregnant women must avoid using this massager.

4. Do not use this massager if you are using

pacemaker or other life support equipment such as an artificial lung or respirator.

5. Should not be used while driving, operating machinery, or during any activity in which involuntary muscle contractions may put the user at undue risk or injury

6. Do not use if you have following medical conditions:

Acute disease, trauma or surgical procedure in past six month

Cardiac dysfunction

Epilepsy

High blood pressure

High fever

cancer

Infectious disease

Malignant tumors

Recent surgery


Skin problems

Pregnancy

Maintenance and care

1. When you are finished using the device, turn it off and take off the pads. Store the unit in a cool, dry place until next use.
2. After usage, please cover the pads with the protective film. Applying a few drops water before and after usage will increase the longevity of the pads.
3. Charge the unit at least every 3 month in event it is not used regularly.
4. Please replace pads after using 15-20 times (Varies widely depending on use, body type and care).
5. Do not apply pads to the same area more than 30 minutes at a time.

Recommend massage points

Troubleshooting

| Problem | Cause | Solution |
|---|---|--|
| During the massage, the skin feels a painful burning sensation or the stimulation becomes weakened. | Adhesive gel pads are not adhering firmly to the skin. Or the gel pads are too dry. | Apply a few drops of water to adhesive surface of each pad and make sure the pads are pressed firmly to the skin during application. |
| Unit is on, but no sensation is felt in the pads. | Both pads are not firmly on the skin. Or there is a loose connection. Or the intensity level needs to be increased. | When you use two pads(one channel) or four pads(two channels) at the same time, make sure you have placed all the pads on your body, otherwise it won't work. Increase the intensity |

| | | |
|--|---|---|
| Adhesive gel pads do not stick to skin even after cleaning and moistening the gel pad. | Adhesive gel pads need to be replaced. | Replace the set of pads. Pads can be purchased from local distributor. |
| Unit does not turn on or the screen is dim. | Battery is low. | Recharge the battery. |
| Does not provide adequate pain relief. | Not using it long enough. Not lying down then massaging neck. | Use your massager 20-30 minutes at a time, 3-6 times per day. Lie on your back for neck support when you need relief for headaches. |

Please Download User Manual on Below Link:

www.moodlightlife.com/downloads

Manufactured by:

Shenzhen Hongqiangxing Electronics Co., Ltd.
Address: 4/F, Jingcheng Building, Xicheng
Industrial District, Xixiang Road, Baoan,
Shenzhen, China

Distributed by:

Wellkang Ltd
ADD: Suite B, 29 Harley Street, London,
W1G9QR, United Kingdom
Support@moodlightlife.com

Made in China

This device complies with part 15 of the
FCC rules

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